

# Colibri

THE COLOURFUL ROUTE TO  
EFFORTLESS, SAFE DRIVING



The Invacare Colibri is a stylish, simple and colourful micro lightweight scooter, designed for those who enjoy a totally independent lifestyle. Less than a minute is all that's needed to disassemble the Colibri for transportation.



Yes, you can.®

**Making Life's  
Experiences Possible**

# Colibri



Thanks to the unique patented Invacare LiteLock™ system the scooter can easily be taken apart without the need for tools and fits neatly into any car boot. Plus, with its small footprint the Invacare Colibri offers greater manoeuvrability in tight, or limited spaces.

Additionally, you can choose from a range of trendy colour panels to match your style, swapping them according to your needs.

Choose where to charge your battery: directly on the Colibri or detached from it.

The soft, padded seat provides comfort during driving. For increased comfort there is also the optional 'easy fit' seat suspension kit.

## Features



















Effortless dismantling can be done in less than a minute for transportation.



A soft padded seat provides additional comfort

## Tech Specs

 <b>Backrest width</b> 465 mm	 <b>Seat depth</b> 400 mm	 <b>Floor to seat height</b> 370 - 420 mm	 <b>Backrest height</b> 385 mm	 <b>Driving unit width</b> 505 mm	 <b>Total length</b> 1010 mm	 <b>Empty weight</b> 49 kg (18 Ah)	 <b>Max user weight</b> 136 kg
 <b>Battery capacity</b> 2 x 18 Ah	 <b>Turning diameter</b> 2200 mm	 <b>Max climable obstacle height</b> 45 mm	 <b>Ground clearance</b> 50mm	 <b>Max safe slope</b> 6° / 10.5°	 <b>Max driving distance</b> Up to 16 km	 <b>Max driving speed</b> 8 km/h	 <b>Motor capacity</b> 200 W

## Colour Options



\* Warranty Terms and Conditions apply. Contact Invacare for full details.

FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email [sales@invacare.co.nz](mailto:sales@invacare.co.nz) or visit [www.invacare.co.nz](http://www.invacare.co.nz)  
**Australia:** Freephone 1800 460 460, email [orders@invacare.com.au](mailto:orders@invacare.com.au) or visit [www.invacare.com.au](http://www.invacare.com.au)  
**Asia:** Phone +61 2 8839 5330, email [anzexportsales@invacare.com](mailto:anzexportsales@invacare.com) or visit [www.invacare.com](http://www.invacare.com)  
**Thailand:** Phone +66 0 2821 5515 or visit [www.invacare.com](http://www.invacare.com)



**Yes, you can.®**